

.

From: Khan, Annie (DPH)
Sent: Tuesday, January 24, 2012 11:21 AM
To: 'laurac@brazoria-county.com'
Subject: Amphetamine Extraction

Hi Laura,

I am not sure how well these extraction methods are going to work with blood samples.

A. Base Conversion (Method 1)

- i. Place about 1-3mg or mL of the sample into a labeled culture tube.
- ii. Add about 1-2mL of 0.1N sodium hydroxide to the sample in the labeled culture tube. Allow the sample to solubilize for about 2-3 hours or overnight.
- iii. Add about 1-2mL of hexane to the labeled culture tube and mix gently.
- iv. Remove the hexane layer from the solution. Transfer the extract into a labeled GC vial and cap tightly.

B. Base Conversion (Method 2)

- i. Place about 1-3mg or mL of the sample into a labeled culture tube.
- ii. Add about 1-4 drops of ammonium hydroxide to the sample in the labeled culture tube.
- iii. Add about 1-2mL petroleum ether to the labeled culture tube and mix gently.
- iv. Remove and transfer the extract into a labeled GC vial and cap tightly.

Hope this helps. Good Luck!!!

Thanks,

Annie
Drug Analysis Lab
617-983-6622 (work)
617-983-6625 (fax)